



aurora
SCIENTIFIC | Performance.
Precision.
Progress.
www.AuroraScientific.com

Parts & Accessories

Aurora Scientific Model: 900B

900C

Software Upgrade for Video Sarcomere Length (VSL) and High-Speed Video Sarcomere Length (HVSL).

901A-CAL

Microscope calibration reticule slide. Used to calibrate the VSL and HVSL software. The spacing of the grating on the reticule is 10 μ m.

901A-TRG

Hardware module to trigger the VSL or HVSL program to begin recording video and Sarcomere Length. Connects through the serial port of the computer.

903A

USB to analog output module that allows sarcomere length values measured by the VSL or HVSL program to be output as a varying voltage signal.