



**aurora**  
SCIENTIFIC | Performance.  
Precision.  
Progress.  
[www.AuroraScientific.com](http://www.AuroraScientific.com)

## **Parts & Accessories**

### **Aurora Scientific Model: 100A**

#### 100-10.6

Replacement 10.6 eV lamp for the Aurora Scientific miniPID. The lamp will ionize the molecules of any odorants with an ionization potential less than 10.6 eV.

#### 100-11.7

Replacement 11.7 eV lamp for the Aurora Scientific miniPID. The lamp will ionize the molecules of any odorants with an ionization potential less than 11.7 eV.

#### 151A

Replacement signal/power cable providing RS-485 serial interface to a 100A digitalPID. The cable is 3 meters in length.

#### 152

Replacement signal/power cable providing RS-232 serial interface to a 100A digitalPID. The cable is 3 meters in length.

#### 200-PUMP

Replacement rotary vane pump for use inside the miniPID control electronics.