



## Parts and Accessories

For Aurora Scientific Model: **1205A**

| Model #                              | 1205A   |
|--------------------------------------|---|
| <b>Apparatus</b>                     |   |
| 800-EA                               | Replacement electrode assembly  |
| 805-LTC                              | Lower tissue clamp assembly for 805A rat <i>in-vitro</i> apparatus – includes: electrode clamp, lower tissue clamp, lower tissue plate, support tube and support tube clamp |
| 800-FT                               | Adaptor mount for 400A series force transducers – converts 800A/805A to strictly isometric  |
| 800-TB-100                           | 100mL Radnoti tissue bath which mounts to 805A rat <i>in-vitro</i> apparatus  |
| 800-OXY                              | Oxygenating tube for Radnoti tissue bath  |
| 800-TC-18                            | Teflon thermocouple probe for measuring bath temperature  |
| <b>Dual-Mode Lever System</b>        |   |
| 300-3B                               | 3 cm long replacement lever arm, attaches to 300C or older 300B series motors for <i>in-vitro</i> and <i>in-situ</i> measurement  |
| 300-4B                               | 4 cm long replacement lever arm, attaches to 305C or older 305B series motors for <i>in-vitro</i> and <i>in-situ</i> measurement  |
| 300-CLA                              | Custom length replacement lever arm, attaches to 300C/305C or older 300B/305B series motors for <i>in-vitro</i> and <i>in-situ</i> measurement                              |
| 300-CBL-3                            | 3m replacement motor cable (Other lengths available on request)   |
| 300-CBL-FR                           | 0.4m long cable to reverse force direction of 300C/305C motor   |
| 6350                                 | Replacement motor for 300C dual-mode lever system (Lever arm sold separately)   |
| 6650-LR                              | Replacement motor for 305C dual-mode lever system (Lever arm sold separately)   |
| <b>Stimulator</b>                    |   |
| 606A                                 | 150cm long BNC to BNC cable   |
| <b>Data Acquisition and Analysis</b> |   |
| 603D                                 | 16-bit NI X Series USB A/D with interface   |
| 603E                                 | 16-bit NI X Series PCI express A/D card   |
| 606A                                 | 150cm long BNC to BNC cable   |
| 607B                                 | 50 Ohm BNC Tee adapter allowing connection of two BNC cables to one BNC port  |
| 615B                                 | 2 <sup>nd</sup> Copy of DMC and DMA suite software  |
| 615C                                 | DMC and DMA suite software upgrade  |
| 608C                                 | 22" LCD monitor   |
| <b>Temperature Control</b>           |   |
| 827A                                 | Heated Water Circulator for maintaining temperature of bath solution (North America power ONLY)   |
| <b>Rack Mounting</b>                 |   |
| 830A                                 | Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking (works with models: 400A series and 701C)   |
| 831A                                 | Heavy gauge, steel desktop or shelf-top rack, with standard 19" width for any Aurora Scientific instruments (6U high for 1200A)   |