



**aurora**  
SCIENTIFIC | Performance.  
Precision.  
Progress.  
[www.AuroraScientific.com](http://www.AuroraScientific.com)

## Parts & Accessories

### **Aurora Scientific Model: 1310A**

808-1

Vertical pin clamp that utilizes a 27G hypodermic needle to immobilize the mouse knee joint during *in-situ* experiments. It mounts on our family of 3-in -1 apparatus.

808-2

Vertical Screw Clamp that uses sharpened, stainless steel screws and a Delrin pad to fix mouse knee joint during *in-situ* experiments. It mounts on our family of 3-in-1 apparatus.

808-3

Upright knee clamp designed for rat that uses a horizontal stainless steel screw to restrain knee joint during *in-vivo* or *in-situ* experiments. It mounts on our family of 3-in-1 apparatus.

808-4

Upright knee clamp designed for rat that uses a dowel pin with a pad on the end for gentle immobilization of rat knee during *in-vivo* or *in-situ* experiments. It mounts on our family of 3-in-1 apparatus.

808-5

Upright knee clamp designed for mouse that uses a horizontal, sharpened stainless steel screw to prevent hind limb/joint movement during *in-vivo* or *in-situ* experiments. It mounts on our family of 3-in-1 apparatus.

808-6

Upright knee clamp designed for mouse that uses a horizontal dowel pin with pad for non-invasive knee restraint during *in-vivo* or *in-situ* experiments. It mounts on our family of 3-in-1 apparatus.

809C-25,

Horizontal muscle bath which replaces the animal platform on 809C apparatus for *in-vitro* intact muscle mechanics experiments. The bath couples to the apparatus by removing the 4 screws for the limb platform and using them to attach the bath in its place. The volume of the bath is 25mL. Tools, tubing and fittings are also provided.

801-RES

1L reservoir assembly complete with built in, oxygenating bubbler. The reservoir is meant to supply oxygenated media to the 809C-25. The reservoir has an outer glass water jacket and can be connected to a circulating water reservoir to maintain perfusate temperature. Tubing and fittings to couple the reservoir are provided.

830A

Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking. Works with models: 200B, 315C/322C, 400A series, 701C and 825A.



**aurora**  
SCIENTIFIC  
[www.AuroraScientific.com](http://www.AuroraScientific.com)

Performance.  
Precision.  
Progress.

## **Parts & Accessories**

### **Aurora Scientific Model: 1310A**

831A

Heavy gauge, steel desktop or shelf-top rack, with standard 19" width for any Aurora Scientific instruments. Available in various heights. (6U high for 1200A and 1300A, 5U high for 1400A, 1500A and 1600A systems).