



**aurora**  
SCIENTIFIC | Performance.  
Precision.  
Progress.  
[www.AuroraScientific.com](http://www.AuroraScientific.com)

## **Parts & Accessories**

### **Aurora Scientific Model: 1600A**

#### **900B**

Video Sarcomere Length (VSL) measurement system includes camera and VSL software that integrates with the 600A Digital Controller. The included USB camera operates at a frame rate of 80 frames per second (FPS) for full frames and up to a maximum of 1000 FPS for reduced frame sizes. The camera mounts via a C-mount adapter (not provided).

#### **901B**

High-Speed Video Sarcomere Length (HVSL) measurement system includes camera and HVSL software that integrates with the 600A Digital Controller. The included Gigabit Ethernet camera operates at a frame rate of 300 frames per second (FPS) for full frames and up to a maximum of 4700 FPS for reduced frame sizes. The camera mounts via a C-mount adapter (not provided).

#### **830A**

Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking. Works with models: 200B, 315C/322C, 400A series, 701C and 825A.

#### **831A**

Heavy gauge, steel desktop or shelf-top rack, with standard 19" width for any Aurora Scientific instruments. Available in various heights. (6U high for 1200A and 1300A, 5U high for 1400A, 1500A and 1600A systems).