



aurora
SCIENTIFIC | Performance.
Precision.
Progress.
www.AuroraScientific.com

Parts & Accessories

Aurora Scientific Model: 300C-I

300-3I

3 cm long replacement lever arm which attaches to the 300C or older 300B series motors. Can be coupled with a short needle or Teflon indentation tips for mechanical stimulation of tissue or biomaterials.

300-4I

4cm long replacement lever arm which attaches to the 300C or older 300B series motors. Can be coupled with a short needle or Teflon indentation tips for mechanical stimulation of tissue or biomaterials.

300-3I-0.2

Teflon tip for use with any indentation arm that is precisely 0.2mm in diameter.

300-3I-0.5

Teflon tip for use with any indentation arm that is precisely 0.5mm in diameter.

300-3I-0.8

Teflon tip for use with any indentation arm that is precisely 0.8mm in diameter.

300-3I-1.0

Teflon tip for use with any indentation arm that is precisely 1.0mm in diameter.

300-3I-1.5

Teflon tip for use with any indentation arm that is precisely 1.5mm in diameter.

300-3I-2.0

Teflon tip for use with any indentation arm that is precisely 2.0mm in diameter.

300-3I-3.0

Teflon tip for use with any indentation arm that is precisely 3.0mm in diameter.

300-CBL-3

Replacement motor cable to connect the motor to the control electronics. The cable is three meters long. (Other lengths available on request).

300-CBL-10

Replacement motor cable to connect the motor to the control electronics. The cable is ten meters long. (Other lengths available on request).

300-CBL-FR

0.4 meter long cable to reverse length and force direction of 300C/305C motors.

6350

Replacement motor for 300C-I, indenter system. (Lever arm sold separately).



aurora
SCIENTIFIC
www.AuroraScientific.com

Performance.
Precision.
Progress.

Parts & Accessories

Aurora Scientific Model: 300C-I

6650-LR

Replacement motor for 305C-I, indenter system. (Lever arm sold separately).