



aurora
SCIENTIFIC
www.AuroraScientific.com

Performance.
Precision.
Progress.

Parts & Accessories

Aurora Scientific Model: 300C

300-FP Footplate for 300C for Mouse
305-FP Footplate for 305C for Rat
305-FP-M Footplate for 305C for Mouse
300-3B Lever Arm for 300C
300-4B Lever Arm for 305C
300-4B-HI Lever Arm for 305C, High Inertia
300-CLA Custom Lever Arm
300-3I Indenter Lever Arm for 300C
300-4I Indenter Lever Arm for 305C

Indenter Tips:

300-3I-0.2, 0.5, 0.8, 1.0, 1.5, 2.0, 3.0

Cables:

300-CBL-3 3m Motor Cable for 300C/305C
300-CBL-10 10m Motor Cable for 300C/305C
300-CBL-FR Force Reversal Cable for 300C/305C

Motors:

6350 for 300C/300C-LR
6650-LR for 305C/305C-LR
6900 for 309C
6400 for 310C/310C-LR

6650-LR for 305C/305C-LR
6900 for 309C
6400 for 310C/310C-LR