



**aurora**  
SCIENTIFIC | Performance.  
Precision.  
Progress.  
[www.AuroraScientific.com](http://www.AuroraScientific.com)

## **Parts & Accessories**

### **Aurora Scientific Model: 315C/322C**

#### **300-1B**

1 cm long replacement lever arm for 315C and older 312B/C high-speed length controllers. Note that the arm must be attached to the slotted shaft with epoxy.

#### **300-2B**

2 cm long replacement lever arm for 322C high-speed length controller and older 318B high speed length controllers. This arm fastens to the shaft using a 1/16" Allen wrench (included).

#### **6215H**

Replacement motor for 315C High Speed Length Controller. (Lever arm sold separately).

#### **6220H**

Replacement motor for 322C High Speed Length Controller. (Lever arm sold separately).

#### **312-CBL-3**

Replacement motor cable to couple the motor to the control electronics. The cable is three meters long. (Other lengths available on request).

#### **830A**

Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking. Works with models: 200B, 315C/322C, 400A series, 701C and 825A.