



aurora
SCIENTIFIC | Performance.
Precision.
Progress.
www.AuroraScientific.com

Parts & Accessories

Aurora Scientific Model: 315C/322C

300-1B

1 cm long replacement lever arm for 315C and older 312B/C high-speed length controllers. Note that the arm must be attached to the slotted shaft with epoxy.

300-2B

2 cm long replacement lever arm for 322C high-speed length controller and older 318B high speed length controllers. This arm fastens to the shaft using a 1/16" Allen wrench (included).

6215H

Replacement motor for 315C High Speed Length Controller. (Lever arm sold separately).

6220H

Replacement motor for 322C High Speed Length Controller. (Lever arm sold separately).

312-CBL-3

Replacement motor cable to connect the motor to the control electronics. The cable is three meters long. Other lengths are available on request.

830A

Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking. Works with models: 200B, 315C/322C, 400A series, 701C and 825A.