



aurora
SCIENTIFIC | Performance.
Precision.
Progress.
www.AuroraScientific.com

Parts & Accessories

Aurora Scientific Model: 400A

400-TR

Repair kit for replacing broken or damaged output tubes. The kit includes: a tube alignment jig, ten output tubes, two-part epoxy, a butane powered heating iron kit, various Allen keys and an X-Acto knife complete with a No. 17 chiselling blade. See our blog for more detailed instructions on use.

400-OT-7

10 pack of glass output tube for replacing broken or damaged tubes on 400A, 403A, 404A, 405A, and 406A transducers models. For the 402A and 407A models the tubes are stainless steel.

830A

Half-rack adapter that allows half-rack instruments to be mounted in standard 19" racking. Works with models: 200B, 315C/322C, 400A series, 701C and 825A.