



Parts and Accessories

For Aurora Scientific Model: 400B

Model #	400B
Controller	
400-AMP	Controller
Heads	
400B-HD	50mN Head
402B-HD	500mN Head
403B-HD	5mN Head
404B-HD	100mN Head
405B-HD	10mN Head
406B-HD	0.5mN Head
407B-HD	1000 mN Head
Output Tube	
400-TR	Repair kit for replacing broken or damaged output tubes that includes a tube alignment jig, ten output tubes, epoxy, butane powered torch, various Allen keys, X-Acto knife and No. 17 chiselling blades (See our blog for more detailed instructions on use)
400-OT-7	10-pack of 7mm glass output tubes for replacing broken or damaged tubes on 400B, 403B, 404B, 405B and 406B models.
400-OT-7-SS	10-pack of 7mm stainless steel output tubes for replacing tubes on 402B and 407B models.
Rack Mounting	
830A	Half-rack adapter that allows 400B half-rack control box to be mounted in standard 19" racking
831A	Heavy gauge, steel desktop or shelf-top rack, with standard 19" width for Aurora Scientific instruments