



aurora
SCIENTIFIC

Performance.
Precision.
Progress.

www.AuroraScientific.com

Specifications

For Aurora Scientific Model: **604A**

Model #	604A
Interface Specifications	
Connector Type	BNC – isolated, male
A/D Inputs	8: Length In, Force In, Aux 1-6 In
D/A Outputs	2: Length Out, Force Out
Digital I/O	2 Inputs: Trg In 1, Trg In 2 4 Outputs: Trg Out 1, Trg Out 2, Stimulator, Inhibit
PC Connection	Ribbon Cable to NI A/D Card (See 603C)
General Specifications	
Weight [kg]	1.7
Dimensions [cm]	42W (19" rack mount) x 13.5D x 4.5H (1U)