



Specifications

For Aurora Scientific Models: 806D

| Model # | 806D |
|-----------------------------------------|----------------------------------------------------------|
| Apparatus | |
| Limb Platform Dimensions [mm] | 82.5 W x 89 L |
| Animal Platform Dimensions [mm] | 190 W x 305 L |
| Platform Horizontal Movement [mm] | 225 |
| Lever System ¹ Movement [mm] | Vertical: 41 Side-to-side: 125 |
| Resting Tension Adjustment [mm] | ±10 |
| Tension Travel per Knob Rotation [mm] | Course: 20 Fine: 2.6 |
| Limb Clamp Fixation | 25 mounting locations |
| Temperature Control | |
| Type | Circulated heated water |
| Fittings | Dry-break, quick connect |
| General Specifications | |
| Weight [kg] | 5.0 |
| Dimensions [cm] | Overall: 31 W x 53 L x 13 H Frame: 28 W x 50 L x 13 H |

¹ Compatible with 300C and 305C lever systems