



## Specifications

**For Aurora Scientific Models: 809C**

Model #	809C
<b>Apparatus</b>	
Animal Platform Dimensions [mm]	90 W x 115 L
Platform Horizontal Movement [mm]	100
Lever System <sup>1</sup> Movement [mm]	Vertical: 41 Side-to-side: 125
Resting Tension Adjustment [mm]	±10
Tension Travel per Knob Rotation [mm]	Course: 20 Fine: 2.6
Limb Clamp Fixation	26 mounting locations
<b>Temperature Control</b>	
Type	Circulated heated water
Fittings	Dry-break, quick connect
<b>General Specifications</b>	
Weight [kg]	1.7
Dimensions [cm]	Overall: 21 W x 28 L x 12 H Frame: 21 W x 26 L x 12 H

<sup>1</sup> Compatible with 300C and 305C lever systems